



Claims

I claim:

Claim 1 (cancelled).

Claim 2 (new). An erect-king band, used for treatment of impotence and premature ejaculation, comprising:

- (a) a rubber band;
- (b) a first hook-and-loop fastener, which is not wider than said band, attached to a first end of said band;
- (c) a second hook-and-loop fastener, which is not wider than said band, attached to a second end of said band and on the surface of said band opposite to the surface which said first fastener is attached to;
- (d) a first marker on said first fastener to indicate the position where the front end of said second fastener reaches when the two said fasteners join together during the first stage of erection;
- (e) a second marker on said first fastener to indicate the position where the front end of said second fastener reaches when the two said fasteners join together during the second stage of erection to achieve full erection status;

wherein two said fasteners join together at the positions defined by two said markers to form a loop structure encircling the penis shaft at the junction of the penis and scrotum wherein the diameter of said loop structure defined by said first marker during the first stage of erection is larger than the one defined by said second marker during the second stage of erection and wherein the elastic force is between 0.35-0.70kg when said loop structures are formed on the penis at the junction of the penis and scrotum.

Claim 3 (new). An erect-king band as defined in claim 2, in which said band further meets

the following dimensions:

- (a) its width is in the range of 12.5-13.5 mm;
- (b) its thickness is in the range of 1.2-1.3 mm;
- (c) its hardness is in the range of 47-53 on the Shore "A" scale;
- (d) its length varies dependent on the diameter of the penis so that it is able to form a loop to encircle the penis when it is extended by 20-30%.

Claim 4 (new). An erect-king band as defined in claim 3, in which said hook-and-loop fasteners further have the following dimensions:

- (a) their thickness is in the range of 1-1.10 mm;
- (b) their width is in the range of 12.5-13.5 mm and not wider than said band;
- (c) their length is in the range of 28-32 mm;
- (d) the distance between said two markers on said first fastener is 1.0-1.7 cm.

Claim 5 (new). A method of using an erect-king band as defined in claim 2 to treat impotence and premature ejaculation, comprising:

- (a) lifting a penis and placing a erect-king band at the junction of the penis and scrotum;
- (b) then stretching said band and encircling the penis;
- (c) and then attaching said second fastener with said first fastener at said first marker;
- (d) and then massaging spermary and scrotum, and pressing the spermary and blood vessel to boost blood to flow passing said band and into the sponge-like tissue until half erection status is achieved;
- (e) and then releasing two said fasteners by pulling them apart;
- (f) and then re-stretching said band and re-encircling the penis;
- (g) and then re-attaching said second with said first fastener at said second marker;
- (h) and then massaging spermary and scrotum, and press the spermary and blood vessel to

boost blood to flow passing said band and into the sponge-like tissue until full erection status is achieved.